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**Bibliographie :** 1. Penninx JP, Mol BW, Bongers MY. Endometrial ablation with paracervical block. J Reprod Med. 2009 Oct;54(10):617-20 2. Penninx JP, Herman MC, Kruitwagen RF, Ter Haar AJ, Mol BW, Bongers MY. Bipolar versus balloon endometrial ablation in the office: a randomized controlled trial. Eur J Obstet Gynecol Reprod Biol. 2016 Jan;196:52-6 3. Reinders IMA, Geomini PMAJ, Leemans JC, et al. Intrauterine fundal anaesthesia during endometrial ablation in the office: A randomised double-blind, non-inferiority trial. Eur J Obstet Gynecol Reprod Biol. 2020 Nov;254:206-211 4. Reinders I, Geomini P, Timmermans A, et al. Local anaesthesia during endometrial ablation: a systematic review. BJOG. 2017 Jan;124(2):190-199 5. Kumar V, Chodankar R, Gupta JK. Endometrial ablation for heavy menstrual bleeding. Womens Health (Lond). 2016 Jan;12(1):45-52 6. Ghoubara A, Gunasekera S, Rao L, Ewies A. Re-intervention and patient satisfaction rates following office radiofrequency endometrial ablation: a comparative retrospective study of 408 cases. J Obstet Gynaecol. 2022 Jul;42(5):1358-1364 7. Ajmi A, Das S. Outpatient endometrial ablation: audit of outcomes and patient satisfaction. J Obstet Gynaecol. 2020 Nov;40(8):1145-1147 8. Kalkat RK, Cartmill RS. NovaSure endometrial ablation under local anaesthesia in an outpatient setting: An observational study. J Obstet Gynaecol. 2011;31(2):152-5 9. Clark TJ, Samuel N, Malick S, Middleton LJ, Daniels J, Gupta JK. Bipolar radiofrequency compared with thermal balloon endometrial ablation in the office: a randomized controlled trial. Obstet Gynecol. 2011 Jan;117(1):109-118. 10. Glasser MH. Practical tips for office hysteroscopy and second-generation "global" endometrial ablation. J Minim Invasive Gynecol. 2009 Jul-Aug;16(4):384-99. 11. Skensved H. Combining Paracervical Block With a Complete Fundal Block Significantly Reduces Patients' Perception of Pain During Radio-Frequency Endometrial Ablation in an Office Setting. J Minim Invasive Gynecol. 2015 Nov-Dec;22(6S):S45 12. Samuel, N.C., Karragiannidou, E. & Clark, T.J. Outpatient versus day-case endometrial ablation using the NovaSure impedance-controlled ablative system. Gynecol Surg 6, 3–9 (2009) 13. Richmond D and Sherwin R. GIRFT Programme National Specialty Report, Maternity and Gynaecology. September 2021. Disponible sur la page <https://www.gettingitrightfirsttime.co.uk/surgical-specialty/obstetrics-and-gynaecology/>. Accès le 21 juin 2022

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MISC-09383-FRA-901 Rév. 001 © 2023 Hologic, Inc. Hologic, NovaSure et les logos associés sont des marques commerciales et/ou des marques déposées d'Hologic, Inc. et/ou de ses filiales aux États-Unis et/ou dans d'autres pays. Le système NovaSure est conçu pour l'ablation de l'endomètre qui recouvre l'intérieur de l'utérus pendant la pré-ménopause chez les femmes souffrant de ménorragie (saignements excessifs) due à des causes bénignes et ne devant plus avoir d'enfant. • Classe IIb • Lire attentivement les instructions figurant dans le manuel d'utilisation. Ces informations sont destinées aux professionnels de santé. Elles ne doivent pas être considérées comme des sollicitations ni faire la promotion de produits lorsque ces activités sont interdites. La documentation Hologic étant distribuée par l'intermédiaire de sites Internet, de diffusions en ligne et de salons professionnels, il n'est pas toujours possible d'en contrôler la disponibilité. Pour obtenir des informations spécifiques sur les produits disponibles à la vente dans un pays donné, contactez votre représentant Hologic local ou écrivez à [france@hologic.com](mailto:france@hologic.com).